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

## Raman Spectroscopy Modular Systems

Flexible and configurable for unique needs

Our modular products allow you to build the exact Raman system you need, while maintaining the flexibility to change and optimize components as your needs evolve.

A modular system begins with your choice of excitation laser, typically routed to the sample with a specialized Raman probe (free-space coupling is also possible). Scattered Raman light is collected by the same probe and routed to a spectrometer that is configured with the sensitivity, optical resolution and Raman shift range needed for your application. We even offer sample holders to facilitate measurement of liquids in cuvettes and vials.

Modular systems offer the most value and flexibility when getting started with a new Raman application, and transition easily into OEM systems. Our modular Raman components include multiple lasers, probes, and user-configurable spectrometers to suit almost any need.

Modular Systems			
Excitation Laser	Multiple options : 532nm, 785nm, custom		
Raman Probe	Multiple options: general purpose, immersion, process		
Spectrometer	QE Pro Series 	Vantana Series 	Maya2000 Pro Series 
Features	* Maximum flexibility * Weak raman signals	* High throughput * Fast measurements	* Cost effective * Broad spectral range
Preconfigured Spectrometers	QE Pro-Raman (785nm)	Ventana-532-Raman Ventana-785-Raman Ventana-785L-Raman	Maya2000 Pro-NIR ( 785nm, 808nm )
Custom Configuration	* User-defined spectral range * Changeable slits to optimize resolution	None	* User-defined spectral range and resolution

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## Modular Spectrometers

Raman spectroscopy requires detection that balances sensitivity and signal to noise with resolution. These spectrometers do so elegantly, each optimized for a particular type of application. Our Maya2000 Pro and QE Pro spectrometers can be custom-configured to your exact excitation wavelength, resolution, range and sensitivity requirements, while their preconfigured versions and our Ventana series are ideal for 532 nm or 785 nm excitation.

### Maya2000 Pro

The Maya2000 Pro is a good general purpose configurable spectrometer for Raman, combining a low stray light optical design with an uncooled back-thinned FFT-CCD detector for enhanced sensitivity that is value-priced. High quantum efficiency in the NIR makes it an excellent choice for long-wavelength Raman.

### QE Pro

This premium configurable spectrometer yields the highest quality Raman spectra, both in resolution and signal to noise. Designed with gold-coated mirrors and a back-thinned FFT-CCD detector cooled to -15 °C, the QE Pro offers a low noise floor even at long integration times. Its ability to deliver sharp peaks from very weak Raman signals makes it an excellent choice when your exact application needs are undefined.

### Ventana Series

The Ventana series of preconfigured spectrometers are designed specifically for Raman spectroscopy. By using a volume phase grating and a lower f/# than our other spectrometers (f/1.3 vs. f/4), they deliver higher throughput for fast measurements. This is further enhanced by control of the detector temperature at +15 °C, keeping the noise floor consistent for applications at short integration times.

	Maya2000 Pro-NIR	QE Pro-Raman	Ventana-785-Raman
Raman Shift Range*	0 - 4200cm <sup>-1</sup>	0 - 2800cm <sup>-1</sup>	200- 2000cm <sup>-1</sup>
Resolution @810nm*	13cm <sup>-1</sup>	11cm <sup>-1</sup>	10cm <sup>-1</sup>
Detector Temperature	Unregulated	-15°C	+15°C
S:N Ratio	450:1	1000:1	550:1
Alternate Configurations	Custom Maya2000 Pro	Custom QE Pro with user-changeable slits	532nm version
Best For	* Broad range * High NIR efficiency * General purpose	* Weak signals * Long integration times * Highest quality data	* Low f/# (high throughput) * Fast measurements

\* Based on 785nm excitation

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





## Integrated Systems

### Turnkey operation for targeted applications

Integrated systems are a convenient solution for applications where your excitation wavelength, Raman shift range, and resolution needs are known and unlikely to change. They operate at the more common laser wavelengths, and use free-space optics in place of fiber optic probes. This can improve sensitivity, and enables unique sampling modes like raster orbital scanning (ROS) or microscopy. Our integrated systems outshine the competition in performance and value, from benchtop systems for the lab to a handheld unit for the field.

		Benchtop		Handheld
Features		* Spatial averaging * Lower average power	* Specific location * High efficiency	* Field use * Matching libraries
Excitation Wavelength (nm)	532	●	●	
	638	●		
	785	●	●	●
	808	●		
System		IDRaman Reader 	IDRaman Micro 	IDRaman Mini 
Applications		* SERS and tagging * Delicate samples * Pharmaceutical * Research and QA labs	* Carbon nanotubes * Semiconductors * Life sciences * Art and forensics	* Explosives, narcotics * Materials inspection * Anti-counterfeiting * Process lines

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## IDRaman Reader

### Benchtop Measurement of Liquids and Solids

The IDRaman reader delivers all the performance of a larger, more expensive system, converting seamlessly between solid and liquid samples. Enabled with raster orbital scanning (ROS), it offers low background, high resolution and sensitive measurements of non-uniform and inhomogeneous samples in a footprint that moves easily around the lab.

An adjustable focus sample holder accommodates vials and cuvettes for study of liquids. The source knob on the laser unit allows you to switch easily to measurement of solid samples, SERS substrates and surfaces in the space below the unit, with a user-exchangeable lens for adjustment of focal length.

- \* Turnkey operation for the lab or process line
- \* Measure from bottom of vial for small volumes
- \* Optional side measurement of cuvettes or vials
- \* Downward facing lens for solid samples

### Choose from 5 High-performance Systems

Systems Available	Raman Shift Range	Resolution
IDRaman reader 532	200 – 3200 $\text{cm}^{-1}$	7 $\text{cm}^{-1}$
IDRaman reader 638	200 – 2000 $\text{cm}^{-1}$	8 $\text{cm}^{-1}$
IDRaman reader 785	200 – 3200 $\text{cm}^{-1}$	10 $\text{cm}^{-1}$
IDRaman reader 785HR	200 – 2000 $\text{cm}^{-1}$	4 $\text{cm}^{-1}$
IDRaman reader 808	200 – 2000 $\text{cm}^{-1}$	4 $\text{cm}^{-1}$

## IDRaman Micro

### Imaging of Specific Locations with High Efficiency

The IDRaman micro delivers spatial resolution, excellent collection efficiency, and high quality spectra at a fraction of the size and cost of a traditional Raman microscope. The superior OneFocus optical system performs Raman measurement and imaging in a single focal plane for optimum sensitivity and ease of use. Created for maximum flexibility, this system is ideal for the research or QA/QC lab. It offers considerably better signal to noise than fiber-coupled systems of comparable cost, and better resolution of weak and closely spaced peaks.

- \* High quality imaging of sample
- \* Lever switches measurement from stage to vial
- \* Measure from bottom of vial for small volumes
- \* Interchangeable objectives
- \* Micro or macro view of a specific location

### Choose from 3 High-performance Systems

Systems Available	Raman Shift Range	Resolution
IDRaman micro 532	200 – 3200 $\text{cm}^{-1}$	7 $\text{cm}^{-1}$
IDRaman micro 785	200 – 3200 $\text{cm}^{-1}$	10 $\text{cm}^{-1}$
IDRaman micro 785HR	200 – 2000 $\text{cm}^{-1}$	4 $\text{cm}^{-1}$

## IDRaman mini

### Handheld Raman System

The IDRaman mini is the smallest, most powerful handheld Raman instrument available today. Designed to make rapid and accurate measurements that identify, authenticate and verify samples, the IDRaman mini is ideal for:

- \* Authentication analysis
- \* Counterfeit detection
- \* Rapid material identification
- \* Verification of incoming materials
- \* Harsh and demanding environments
- \* In-line or at-line testing

From rugged field measurements of chemical and explosive agents to quality assurance and quality control in the laboratory, the IDRaman mini is a truly compact choice for fast and accurate measurements.

Engineering Specifications	IDRaman mini
<b>Physical</b>	
Dimensions: (LWH) cm and inches:	9.4 cm x 7.4 cm x 3.8 cm (3.7" x 2.9" x 1.5")
Weight:	0.33 kg (11.75 oz.)
<b>System</b>	
Raman laser excitation wavelength:	785 nm +/- 0.5 nm, 2 cm <sup>-1</sup> line width, stability <0.1 cm <sup>-1</sup>
Laser output power:	100 mW (50 mW at sample)
Range:	400-2300 cm <sup>-1</sup>
Resolution:	12-14 cm <sup>-1</sup>
Acquisition speed:	Visual confirmation in <9.0 seconds
Signal to Noise:	1000:1
Raster Orbital Scanning:	Sampling function allows detection of inhomogeneous samples without sample degradation
Detector:	2048-element back-thinned, NIR -enhanced CCD array
Collection optics:	NA = 0.50; 8 mm working distance and 2.5 mm interrogation area
Sampling method:	vial (2 mL) or point and shoot accessory for solid sampling
Operating temperature:	-20°C to +40°C
Display:	2.8" resistive touch display
Battery life:	>11 hours with 2 AA batteries
Raman library:	Basic library of ~30 compounds (included)
Accessories:	Box of 144 x 2 mL vials (included); point-and-shoot accessory for solid sampling

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## ACCUMAN PRO-500

### Portable Raman Spectrometer

The pharmaceutical industry constantly needs to meet ever-growing compliance and quality standards. Achieving this often brings substantial costs. Many processes require identification of Active Pharmaceutical Ingredients (API) and excipients, often involving sophisticated sample preparation, lengthy testing time, and substantial labor. Opening up packaging also introduces the risk of contamination and exposure. Both risks and costs could be greatly reduced by using a non-destructive, point-and-shoot, instantaneous identification solution, directly at the production site.

The ACCUMAN PR-500 is a portable Raman spectrometer, offering rapid, through-packing, accurate identification of pharmaceutical materials. The Raman spectrometer quickly captures the molecular fingerprint of target materials, enabling packet-by-packet or batch-by-batch verification, while delivering the same accuracy as the conventional procedures. At the core of the ACCUMAN PR-500 is the QE-Series spectrometer, which for its superior performance and quality, was used in NASA's projects. This revolutionary solution brings the answer to the growing compliance requirements in the pharmaceutical manufacturing industry.

### Advantages

**\* Speedy**

High stability and performance integration system, delivering accurate results in seconds, significantly raising efficiency.

**\* Accurate and Reliable**

The core is the Ocean Optics QE-Series spectrometer. Its industry-leading signal-noise ratio and stability allows accurate capturing of molecular "fingerprint", giving highly reliable identification results.

**\* User-friendly**

An ergonomically-designed hand-held probe weights just 330g, and comes with a high-resolution, touchscreen interface, allowing intuitive, single-hand operation.

**\* Customized Library**


Library can be customized according to users' raw material from different vendors, reducing variations and false from standard libraries.

**\* Compliance**

Designed according to GAMP 5, complies with CFR part 11 and GxPs.

User Interface

PC User Interface

### Specification

Spectral range	200~2000 $\text{cm}^{-1}$	200~3000 $\text{cm}^{-1}$	200~3900 $\text{cm}^{-1}$
Wavenumber resolution	4~6 $\text{cm}^{-1}$	5~9 $\text{cm}^{-1}$	10~12 $\text{cm}^{-1}$
Laser wavelength	785nm $\pm$ 0.5nm		
Laser power	350 mW Max, Adjustable in 10 steps		
Optical parameter	NA 0.22, 7.5mm working distance		
Detector	TE cooling Backthinned CCD array		
Signal-Noise-Ratio	10000:1		
Intergration time	0.1~30s		
Screen	4.5"LCD 720p Touch screen		
Bar code	1D / 2D		
Data transfer	USB 2.0		
Data format	.pdf / .csv / .txt		
Battery	Li-ion Battery , duration time > 5h		
Power adapter	100~240V AC 50/60 Hz		
Weight	(Main unit) 3.8kg , (Probe) 330g		
Dimensions	(Main unit) 29cmx22cmx10cm , (Probe) 15.5cmx7.4cmx2.5cm		
Operating temperature	0~45°C		
Accessory	Probe cap , calibration cap		
Compliance	USP , EP		
Designed according to GAMP 5 guidance , complies with CFR Part 11 and GxPs computerized system requirements.			

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
## QE Pro-Raman

### High-Sensitivity Spectrometer for Raman

#### Wide Dynamic Range

The QE Pro-Raman is a preconfigured for 785 nm Raman excitation, using a 900 lines/mm grating and 50  $\mu\text{m}$  slit. To complete your modular Raman system, add a 785 nm Raman excitation laser, sample holder, software and other Raman accessories.

The high-sensitivity QE Pro is a good option for measuring the Raman shift of aromatic hydrocarbons such as toluene.

Engineering Specifications	QE Pro-Raman
<b>Physical</b>	
Dimensions: (L x W x H) mm and inches	182 x 110 x 47 mm (7.17" x 4.33" x 1.85")
Weight: kg and lb	Spectrometer: 1.15 kg (2.6 lbs.); power supply: 0.45 kg (1 lb.)
<b>Detector</b>	
Type:	Hamamatsu scientific grade, back-thinned, TE Cooled, 1044 x 64 element CCD array
Range:	185-1100 nm
Quantum efficiency:	90% (peak)
<b>Spectroscopic</b>	
Wavelength range:	780-1000 nm
Raman shift:	0-2800 $\text{cm}^{-1}$
Resolution:	7 – 11 $\text{cm}^{-1}$ over the range*
Integration time:	8 ms to 60 minutes
Dynamic range:	~85,000:1
Signal-to-noise ratio:	System: 1000:1 (single acquisition)
Stray light:	<0.08% at 600 nm; 0.4% at 435 nm
Buffering:	15,000 spectra
Fiber optic connector:	SMA 905 and Ocean Optics FC
<b>Electronics</b>	
Power requirement:	Supply voltage: 4.5 – 5.5 V
Strobe functions:	continuous and single strobe
Interfaces:	USB 2.0, 480 Mbps (USB 1.1 compatible); RS-232 (5-wire)
Temperature:	TE Cooler can only cool 40 $^{\circ}\text{C}$ below ambient temperature; Operation: -40 $^{\circ}\text{C}$ to +50 $^{\circ}\text{C}$
Humidity:	$\leq$ 90% noncondensing

#### Example Configurations and Predicted Performance

Grating	Slit	Raman Shift	Resolution
1200 lines/mm	10 $\mu\text{m}$	150-4000 $\text{cm}^{-1}$	~8 $\text{cm}^{-1}$
1200 lines/mm	25 $\mu\text{m}$	150-4000 $\text{cm}^{-1}$	~10 $\text{cm}^{-1}$
600 lines/mm	10 $\mu\text{m}$	150-7500 $\text{cm}^{-1}$	~16 $\text{cm}^{-1}$
600 lines/mm	25 $\mu\text{m}$	150-7500 $\text{cm}^{-1}$	~19 $\text{cm}^{-1}$

#### QE Pro with 785 nm Raman Laser Excitation

Grating	Slit	Raman Shift	Resolution
1200 lines/mm	50 $\mu\text{m}$	150-2100 $\text{cm}^{-1}$	~6 $\text{cm}^{-1}$
1200 lines/mm	100 $\mu\text{m}$	150-2100 $\text{cm}^{-1}$	~8 $\text{cm}^{-1}$
600 lines/mm	50 $\mu\text{m}$	150-3950 $\text{cm}^{-1}$	~13 $\text{cm}^{-1}$
600 lines/mm	100 $\mu\text{m}$	150-3950 $\text{cm}^{-1}$	~18 $\text{cm}^{-1}$

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## VENTANA Series Spectrometer

### High-Sensitivity Spectrometer for Raman Measurements


VENTANA spectrometers are high throughput, high sensitivity benchtop instruments for your lab and QA/QC environments. They are a great choice for Raman analysis and fluorescence measurements, offering short integration times for fast measurements.

Specifications	VENTANA 532 Raman	VENTANA 785 Raman	VENTANA 785L Raman	VENTANA VIS-NIR
				
<b>Physical</b>				
Dimensions:	170 x 110 x 50 mm	125 x 110 x 50 mm	127 x 159 x 51 mm	170 x 110 x 50 mm
Weight:	0.9 kg (2 lb.)			
<b>System</b>				
Wavelength range:	533-690 nm	800-940 nm	800-940 nm	430-1100 nm
Raman shift range:	350-4300 $\text{cm}^{-1}$	200-2000 $\text{cm}^{-1}$	200-2000 $\text{cm}^{-1}$	None
Resolution (FWHM):	20 $\text{cm}^{-1}$	10 $\text{cm}^{-1}$ @ 810 nm	10 $\text{cm}^{-1}$ @ 810 nm	4.0 nm
Signal to Noise (typical):	~550:1			
Dynamic range:	17000:1			
Thermoelectric cooling:	None	15 °C	15 °C	None
Operating environment:	0-50 °C (non-condensing)	0-50 °C (noncondensing)	5-45 °C (noncondensing)	0-50 °C (non-condensing)
Fiber optic input:	Optimized for 600 $\mu\text{m}$ , NA 0.39 NA; SMA 905			
Power input:	12 VDC			
Interface:	USB 2.0 (no external accessory connector access)			
Operating software:	OceanView (OmniDriver and SeaBreeze device drivers)			
Regulatory:	CE Mark, RoHS			

### Order Information

Model	Description	Example Applications
VENTANA-532-RAMAN	High sensitivity, preconfigured for 532 nm excitation Raman measurements	* Life sciences * Clinical diagnostics
VENTANA-785-RAMAN	High sensitivity, preconfigured for 785 nm excitation Raman measurements	* Pharmaceutical R&D and QA/QC
VENTANA-785L-RAMAN	High sensitivity, preconfigured for 785 nm excitation Raman measurements with integrated 785 nm laser	* Free-space Raman analysis of powders and more
VENTANA-VIS-NIR	High sensitivity, preconfigured with 430 - 1100 nm bandwidth, 50 micron slit	* Life sciences

## Apex 785 Raman Spectrometer

The Apex 785 Raman spectrometer was designed for maximum resolution and sensitivity and delivers exceptional performance. Using patented high throughput virtual slit (HTVS) technology licensed from Tornado Spectral Systems, Apex eliminates the trade-offs among throughput, resolution and spectral range and provides high sensitivity, with excellent optical resolution, over a wide Raman shift range.

HTVS technology dramatically increases the throughput of the detection system, allowing faster measurements and lower excitation powers for delicate samples. Measure samples with very low noise in a fraction of the integration time.

Apex is the first of many new spectrometers, light sources and sampling accessories in our Elite series. Elite products deliver superior performance in a compact and customizable format. Achieve the highest quality data, choose Apex.

### Features

- \* High throughput virtual slit for shorter integration times, faster measurements
- \* High resolution and high sensitivity eliminates typical performance trade-offs
- \* Flexible use with a variety of lasers and accessories

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

## NeoFox Phase Measurement System

Benchtop and Handheld Electronics for Your Optical O2 Sensor

For oxygen sensing, the NeoFox® Phase Measurement System is our most popular fluorescence-based optical sensing system. Because of NeoFox's unique ability to improve overall system stability and make calibration easier for a wide variety of oxygen sensing setups, it is the choice for measuring fluorescence lifetime, phase and intensity. Plus, NeoFox is brilliantly suited for applications where sensitivity to drift and system stability are critical.

The NeoFox uses LED excitation and photodiode detection with filter-based wavelength selection for easy experimental setup and control. Because the unit is self-contained, it is invariant to fiber bending and stray light, and has a wide dynamic range of optical intensity as well as low optical and electronic crosstalk, and low drift and phase noise.

NEOFOX-GT is a more robust benchtop option with RS-232 communications capability.

### NeoFox Specifications

	Benchtop Item Code: NEOFOX
Dimensions	107.95 x 63.5 x 38.1 mm
Weight	642 g
Principle	Photoluminescence quenching using a ruthenium compound; sensor measures O2 partial pressure
Parameters measured	Luminescence phase shift, AC luminescence intensity, temperature (via optional external thermistors) and pressure (via onboard pressure transducer)
Sensor coating formulations (sol gel-embedded dyes)	General purpose (FOXY), high-sensitivity (FOSPOR) and hydrocarbon-ready (HIOXY)
Media	Gases and liquids
Computer interface	PC
Operating systems	Windows XP (32-bit); also, Windows 7 (32-bit)
Power input	5VDC, 500 mA steady state
Communications	USB, analog out, RS-232

### NeoFox Viewer Software

NeoFox Viewer is the Windows-based software that allows you to collect, manage and analyze data with your NeoFox measurement system. This nimble software also makes it simple to configure your NeoFox and update firmware when necessary. Use NeoFox Viewer to get the most out of your NeoFox and to ensure reliability through each measurement.

## NeoFox Probe and Patch-Based Systems

### Full System Performance Specifications

One of the biggest advantages our fiber optic oxygen sensors offer compared with electrodes and other optical sensors is the range of available sampling options. In the table below we've compiled performance data for NeoFox-based systems used with our two primary sampling choices: oxygen probes and patches. Please note that some performance parameters vary according to the sensor format and coating formulation used.

One other item of interest: The specifications for O<sub>2</sub>% and dissolved oxygen range are given for conditions at 1ATM (atmosphere). At 1 ATM (typical conditions on Earth), we breathe many molecules such as nitrogen, helium, hydrogen and oxygen. All these molecules make up the total pressure in the environment. The ruthenium in our sensors is sensitive only to oxygen, which is just one part of the total pressure in the environment.

So, when we determine the percentage of oxygen present in the environment at 1 ATM, we're measuring part of the total pressure, hence the term "partial pressure."

Probe-based System Specifications	FOXY Formulation	FOSPOR Formulation	HIOXY Formulation
Recommended use	General purpose coating	High-sensitivity coating for lowoxygen environments	Robust coating for hydrocarbon-rich environments
O <sub>2</sub> % range (at 1 ATM)	0-100%	0-10%	0-20.9%
DO range (ppm at 1 ATM)	0-40 ppm	0-4 ppm	0-8 ppm
Temperature range	-50-+80 °C for probes	0-+60 °C for probes	-50-+60 °C for probes
O <sub>2</sub> % resolution	100 ppm in gas	10 ppm in gas	100 ppm in gas
DO resolution (at room temp)	4 ppb	0.4 ppb	4 ppb
O <sub>2</sub> % accuracy	<5% of reading	<5% of reading	<5% of reading
DO accuracy	<5% of reading	<5% of reading	<5% of reading
Min. detectable level in gas	0.01% - 0.05%	0.001% - 0.01%	0.01% - 0.05%
Response time	<1 s in gas	<30-60 s in gas	<1 s in gas
	45-60 s with overcoating in gas	60-90 s with overcoating in gas	NA
	30-45 s in pure water	60-90 s in pure water	~45 s in pure water
Patch-based System Specifications	FOXY Formulation	FOSPOR Formulation	HIOXY Formulation
Recommended use	General purpose coating	High-sensitivity coating for lowoxygen environments	Robust coating for hydrocarbon-rich environments
O <sub>2</sub> % range (at 1 ATM)	0-100%	0-10%	0-20%
DO range (ppm at 1 ATM)	0-40 ppm	0-4 ppm	0-8 ppm
Temperature range	-20 to +60 °C for patches	0 to +60 °C for patches	-20 to +60 °C for patches
O <sub>2</sub> % resolution	0.05% (at 20 s averaging)	0.01% (at 30 s averaging)	0.05% (at 20 s averaging)
DO resolution (at room temp)	20 ppb	4 ppb	20 ppb
O <sub>2</sub> % accuracy	5% of reading	5% of reading	5% of reading
DO accuracy	5% of reading	5% of reading	5% of reading
Min. detectable level	0.1% O <sub>2</sub>	0.01% O <sub>2</sub> (at 30 s averaging)	0.1% O <sub>2</sub>
Min. detectable level in water (at room temp)	40 ppb	4 ppb	40 ppb
Response time	<1 s in gas	30-60 s	<1 s in gas
	~30-45 s with overcoating in gas	~60-90 s with overcoating in gas	NA
	~45 s in pure water	~60-90 s in pure water	~30-45 s in pure water
Stability (Continuous LED)	FOXY Formulation	FOSPOR Formulation	HIOXY Formulation
Lifetime stability (Tau)	0.0006 μsec/hour	0.003 μsec/hour	0.0002 μsec/hour
Oxygen stability %	0.01% hour	0.005% hour	0.007% hour
Modulation range	0.73 kHz-93.75 kHz	0.73 kHz-93.75 kHz	0.73 kHz-93.75 kHz

## FoxyKits for Oxygen Sensing

Fully Integrated Systems for Your Probe or Patch Applications

We've packaged everything you'll need for probe- or patch-based oxygen sensing applications into two convenient kits.

NEOFOX-KIT-PROBE and NEOFOX-KIT-PATCH are complete, out-of-the-box solutions for a variety of benchtop applications in research environments, teaching labs and commercial labs. Each kit allows you to set single- or multi-point calibration and to display oxygen readings in percent oxygen, partial pressure, moles per liter and more.

You'll note that we selected FOXY sensor formulation for both the NeoFox probe and patch kits. FOXY is our most versatile sensor coating chemistry and is useful for a wide range of applications. If FOXY is not appropriate for your applications, you can always opt for one of our modular oxygen sensing systems with a different coating formula.

Here's what you get with each kit:

### NEOFOX-KIT-PROBE

- \* NeoFox Phase Fluorometer
- \* NeoFox Viewer software
- \* NeoFox-TP temperature probe (thermistor)
- \* 1000  $\mu\text{m}$  bifurcated optical fiber and SMA 905 splice bushing for coupling NeoFox to the probe
- \* FOXY-R multipurpose probe

Probe kits are great for applications involving biological samples such as tissue and organic matter, foods and beverages and liquids in natural environments.

### NEOFOX-KIT-PATCH


- \* NeoFox Phase Fluorometer
- \* NeoFox Viewer software
- \* NeoFox-TP temperature probe (thermistor)
- \* 1000  $\mu\text{m}$  bifurcated borosilicate optical fiber
- \* 5-unit package of 8 mm diameter, self-adhesive FOXY sensor patches

Patch kits are useful anywhere non-intrusive and through-package oxygen measurements are necessary: headspace in food packaging, process monitoring in bioreactor environments and partial pressure of oxygen in biomedical vessels.

## Oxygen Sensor Probes

### Options for a Wide Range of Sensing Environments

Oxygen sensor probes typically are available in a variety of designs and with each of our standard coating formulations (FOXY, FOSPOR and HIOXY). Custom probes and accessories are also available. In most cases, you'll need a 21-02 Splice Bushing and a BIF-BORO bifurcated optical fiber to couple your probe to the NeoFox electronics. Information on the proper care of probes is available elsewhere in this section.

### HypoTube Probes

HypoTube probes are ideal for penetrating septa in packaging, vials and other vessels. The probes are suitable for use in solutions and headspace.

Item	HypoTube Oxygen Probe
Use	Puncturing septa, packaging and other rigid materi- als; especially good for solutions
Core diameter	300 μm
Outer diameter	1000 μm
Length	2 m (assembly), 127 mm (tube)
Ferrule / jacketing	PVC Monocoil
Reconditioning available	Yes
Options	FOXY-HPT-1-PNA general-purpose FOSPOR-HPT-1-PNA high-sensitivity HIOXY-HPT-1-PNA hydrocarbon-ready HPT-PNA (puncture accessory only)

### Electrode-replacement Probes

Our OR125-series probes are designed to replace standard 1/8" (0.125 mm) electrode probes. We offer smooth and O-ring grooved versions.

Item	Direct-replacement Probes for O-ring grooved electrodes
Use	Replacement for 1/8" OD (3.175 mm) oxygen electrodes (O-ring grooved)
Core diameter	1000 μm
Outer diameter	3.175 mm
Length	63.5 mm
Ferrule / jacketing	Stainless steel, titanium and PEEK versions; smooth or O-ring grooved
Reconditioning available	Yes
Options	FOXY-OR125-G general-purpose FOSPOR-OR125-G high-sensitivity HIOXY-OR125-G hydrocarbon-ready FOXY-OR125-GT general-purpose FOSPOR-OR125-GT high-sensitivity HIOXY-OR125-GT hydrocarbon-ready

Item	Direct-replacement Probes for 1/8" OD (3.175 mm) electrodes
Use	Replacement for 1/8" OD (3.175 mm) oxygen elec- trodes (smooth)
Core diameter	1000 μm
Outer diameter	3.175 mm
Length	63.5 mm
Ferrule / jacketing	Stainless steel, titanium and PEEK versions; smooth or O-ring grooved
Reconditioning available	Yes
Options	FOXY-OR125 general-purpose FOSPOR-OR125 high-sensitivity HIOXY-OR125 hydrocarbon-ready

## Oxygen Sensor Probes

### Options for a Wide Range of Sensing Environments

#### Polyimide Probes

Choose a polyimide probe for applications where the sample environment is hostile to metallic probes. Polyimide also offers good resistance to harsh chemicals. Note: The 200  $\mu\text{m}$  version of these probes is less robust and should be handled carefully. Contact us about custom probe lengths and availability of HIOXY formula probe coatings.

Item	Polyimide Probes
Use	Environments hostile to metallic probes
Core diameter	200 $\mu\text{m}$ and 600 $\mu\text{m}$ versions available
Outer diameter	710 $\mu\text{m}$
Length	2 m (other lengths available)
Ferrule / jacketing	Silicone
Reconditioning available	Yes
Options	FOXY-PI600 and FOXY-PI200 general-purpose FOSPOR-PI600 and FOSPOR-PI200 high-sensitivity

#### General-purpose Probes

Stainless steel 1/16" (1.587 mm) OD probes are versatile options for a range of lab and other applications.

Item	General-purpose 1.587 mm (1/16") Probes
Use	General purpose
Core diameter	1000 $\mu\text{m}$
Outer diameter	1.587 mm (1/16")
Length	152.4 mm
Ferrule / jacketing	Stainless steel
Reconditioning available	Yes
Options	FOXY-R general-purpose FOSPOR-R high-sensitivity HIOXY-R hydrocarbon-ready

Item	General-purpose 1.587 mm (1/16") Probes; short-length versions
Use	General purpose
Core diameter	600 $\mu\text{m}$
Outer diameter	1.587 mm (1/16")
Length	32 mm (other lengths available)
Ferrule / jacketing	Titanium
Reconditioning available	Yes
Options	FOXY-600-32MM general-purpose FOSPOR-600-32MM high-sensitivity HIOXY-600-32MM hydrocarbon-ready

#### Process-ready Probes

Robust 1/4" (6.35 mm) stainless steel probes have a high pressure rating for process environments

Item	Process-ready 6.35 mm probes
Use	Process environments
Core diameter	1000 $\mu\text{m}$
Outer diameter	6.35 mm
Length	177.8 mm
Ferrule / jacketing	Stainless steel
Reconditioning available	Yes
Options	FOXY-T1000 general-purpose FOSPOR-T1000 high-sensitivity HIOXY-T1000 hydrocarbon-ready

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## Oxygen Sensor Probes

### Options for a Wide Range of Sensing Environments

#### Small-Diameter Probes

Slender, aluminum-jacketed probes work well where sampling space is limited. Standard and tissue-monitoring versions are available.

Item	Small-diameter, Al-jacketed Probes
Use	Fine spatial resolution applications
Core diameter	300 $\mu\text{m}$
Outer diameter	500 $\mu\text{m}$
Length	1 m (other lengths available)
Ferrule / jacketing	Aluminum
Reconditioning available	Yes
Options	FOXY-AL300 general-purpose FOSPOR-AL300 high-sensitivity HIOXY-AL300 hydrocarbon-ready

Item	Small-diameter Probe for Tissue Monitoring
Use	Tissue monitoring
Core diameter	300 $\mu\text{m}$ (fiber)
Outer diameter	500 $\mu\text{m}$ (fiber) 25.4 mm (1") at tip, body 1.6 mm (1/16") stainless steel
Length	1 m (other lengths available)
Ferrule / jacketing	Aluminum
Reconditioning available	No
Options	FOXY-AL300-TM general-purpose

#### Respiration Probe

This plastic probe is designed for monitoring oxygen respiration.

Item	Respiration Probes (Standard)
Use	Monitoring of oxygen tension in respiratory gases
Core diameter	200 $\mu\text{m}$
Outer diameter	6.35 mm
Length	107.9 mm
Ferrule / jacketing	Plastic
Reconditioning available	No (uses replaceable glass fiber membranes)
Options	FOXY-RESP general-purpose FOSPOR-RESP high-sensitivity

## Oxygen Sensor Temperature Compensation

### Options for Dealing with Temperature Effects in O2 Measurements

Our optical O<sub>2</sub> sensors are affected by temperature. Temperature affects both fluorescence intensity and excited state lifetime due to the decrease of fluorescence energy quantum efficiency with temperature increase. The effect of change in temperature is seen as a change in the calibration slope.

As temperature increases and partial pressure of oxygen remains the same, compensating for temperature mitigates the risk of a false partial pressure reading. For best results, the sample must be held at a constant temperature ( $\pm 1^\circ\text{C}$ ). The temperature response of our sensors can be determined by the user, or can be supplied by a factory calibration. In a multichannel sensor setup, the user must account for temperature effects in each channel.

### In-house Temperature Calibration Services

If your sample cannot be maintained at a constant temperature ( $\pm 1^\circ\text{C}$ ), you can perform a temperature calibration in NeoFox Viewer software or include our temperature calibration service as part of your order. The standard service covers environments from 0-80°C; extended-range service is available for environments outside those parameters. You'll need to tell us which oxygen sensor formulation you're using (FOXY, FOSPOR or HIOXY) and provide the temperature and oxygen concentration ranges of your sample environment.

Item codes: xxxx (sensor formulation)-CAL (standard option) and xxxx (sensor formulation)-CAL-EXT (extended temperature range)

### Temperature Probe Options for Optical Oxygen Sensor Systems

We offer temperature probe options for use with our oxygen sensor systems: a rugged, discrete thermistor or a T1000-style 1/4" (6.35 mm) stainless steel probe with an embedded thermistor. The latter is available in each of three coating formulations (FOXY, FOSPOR, HIOXY).

T1000-TS-6CM-NEO

T1000-TS-NEO

NEOFOX-TP

### Temperature Probes for Oxygen Sensor Systems

Item	Description	Core Diameter	Outer Diameter	Length	Ferrule/Jacketing	Pressure Rating	Reconditioning Available
FOXY-T1000-TS-NEO	1/4" process probe w/general-purpose coating and thermistor	1000 $\mu\text{m}$	6.35 mm	177.8 mm	Stainless steel	3000 psi	Yes
FOXY-T1000-TS-6CM-NEO	1/4" process probe w/general-purpose coating and thermistor	1000 $\mu\text{m}$	6.35 mm	60 mm	Stainless steel	3000 psi	Yes
FOSPOR-T1000-TS-NEO	1/4" process probe w/high-sensitivity coating and thermistor	1000 $\mu\text{m}$	6.35 mm	177.8 mm	Stainless steel	3000 psi	Yes
FOSPOR-T1000-TS-6CM-NEO	1/4" process probe w/high-sensitivity coating and thermistor	1000 $\mu\text{m}$	6.35 mm	60 mm	Stainless steel	3000 psi	Yes
HIOXY-T1000-TS-NEO	1/4" process probe w/hydrocarbon-ready coating and thermistor	1000 $\mu\text{m}$	6.35 mm	177.8 mm	Stainless steel	3000 psi	Yes
HIOXY-T1000-TS-6CM-NEO	1/4" process probe w/hydrocarbon-ready coating and thermistor	1000 $\mu\text{m}$	6.35 mm	60 mm	Stainless steel	3000 psi	Yes
NEOFOX-TP	3" x #32 AWG thermistor w/10,000 Ohm resistance @ 25 $^\circ\text{C}$	115 mm	3.5 mm	76.2 mm	Radial leaded wires	Not tested	No

## RedEye Optical Sensing Patches

### Self-adhesive Patches for Non-intrusive Oxygen Measurements

The RedEye® indicator patch measures oxygen non-invasively in sealed packaging and containers used in medical, pharmaceutical and food applications. Using a combination of proprietary sensing material and measuring technologies, this non-invasive patch enables quick determination of the presence of oxygen, as well as quantitative measurements. Oxygen monitoring can also ensure patient safety in point-of-care analysis or indicate a sterile seal on surgical instruments and drug packaging. RedEye patches meet USP Class VI certification for biocompatibility.

### RedEye Features a Proprietary Sol Gel Coating

RedEye patches are unique in that high-performance sol gel coatings are used – rather than polymer membranes. RedEye coatings are capable of monitoring low levels of oxygen in gas (to 0.01%) and dissolved oxygen in liquids (to 4 ppb), as well as the higher oxygen levels present in cell culture and respiratory monitoring. Also available is a new design that uses the HIOXY formulation on a BK7 glass substrate autoclavable at 100 °C. Contact an Applications Scientist for details on item RE-HIOXY-HTC.

The RedEye can be integrated into packaging for continuous monitoring or used externally for post-production and R&D monitoring purposes. Depending on the application, the simple presence of oxygen can be visually determined by color change with a handheld LED. A fluorometer can also be used to directly measure oxygen partial pressure.

### Typical Applications

- \* Point-of-care analysis (e.g., disposable oxygen attachments for ventilators used during anesthesia operation)
- \* Blood bag analysis
- \* Beverage and food packaging
- \* Bioprocess control
- \* Cell culture monitoring

Specifications	FOXY Formulation	FOSPOR Formulation	HIOXY Formulation	HIOXY High-Temp Option
Recommended use	General purpose coating	High-sensitivity coating for low-oxygen environments	Robust coating for hydrocarbon-rich environments	Applications where autoclaving is necessary
O2% range (at 1 ATM)	0 - 100%	0 - 10%	0 - 20%	0 - 20.9%
DO range (ppm at 1 ATM)	0 - 40 ppm	0 - 4 ppm	0 - 8 ppm	0 - 8 ppm
Temperature range	0 to +60 °C for patches	0 to +60 °C for patches	0 to +60 °C for patches	0 to +100 °C for glass window (patch)
O2% resolution	0.05% (at 20 s averaging)	0.01% (at 30 s averaging)	0.05% (at 20 s averaging)	0.05% (at 20 s averaging)
DO resolution (at room temp)	20 ppb	4 ppb	20 ppb	4 ppb
O2% accuracy	5% of reading	5% of reading	5% of reading	<5% of reading
DO accuracy	5% of reading	5% of reading	5% of reading	<5% of reading
Min. detectable level in gas	0.1% O2	0.01% O2 (at 30 s averaging)	0.1% O2	0.01% - 0.05%
Min. detectable level in water (at room temp)	40 ppb	4 ppb	40 ppb	40 ppb
Response time	<1 s in gas	30 - 60 s in gas	<1 s in gas	5 - 10 s in gas
	~30 - 45 s with overcoating in gas	~60 - 90 s with overcoating in gas	~30 s with overcoating in gas	~30 s with overcoating in gas
	~45 s in pure water	~60 - 90 s in pure water	~30 - 45 s in pure water	~40 - 60 s in pure water
Patch material	Acrylate	Acrylate	Acrylate	BK7 glass
Patch dimensions	4 mm, 8 mm and 127 mm disk (standard); custom sizes also available	4 mm, 8 mm and 127 mm disk (standard); custom sizes also available	4 mm, 8 mm and 127 mm disk (standard); custom sizes also available	4 mm (standard); custom sizes also available
Standard patch options	Single patch or pack of 5	Single patch or pack of 5	Single patch or pack of 5	Single glass window
Patch-cuvette options	Patch applied to quartz or polystyrene cuvettes	Yes	Yes	No
Overcoat option	Medical-grade overcoat	Medical-grade overcoat	Medical-grade overcoat	Medical-grade overcoat
Adhesive pH compatibility	Yes (pH 4.0 - 10.0)	Yes (pH 4.0 - 10.0)	Yes (pH 4.0 - 10.0)	Yes (pH 4.0 - 10.0)

## TRANS-PH-KIT Series

### Transmissive pH Dip Probe Kits

Engineering Specifications	PH-BCG-TRANS
pH Range	5 – 9
Analytical wavelength	620 nm
Baseline correction wavelengths	509 nm (isosbestic), or > 750 nm
Accuracy	0.05 pH unit (<1% of reading across range)
Resolution	0.02 pH unit
Response time	10 seconds
Calibration requirements (minimum)	Calibration reset: 3 buffers Full calibration: 6 buffers
Sterilization	Gamma radiation, EtO (ethylene oxide)
Chemical compatibility*	Aqueous, alcohols, some organic solvents, peroxides, ammonia, sodium hypochlorite
Discrete stability(lifetime)	50 uses or more. Dispose when absorbance at pH 11 < 0.1 (assumes pH 1 reference)
Storage conditions	Dry storage preferred

### Stainless Steel pH Probe

Engineering Specifications	TRANS-PH-KIT-SS
Probe diameter	6.35 mm (1/4")
Probe length	127 mm
Fiber length	2 meters
Breakout	1.5 meters from the end of the probe
Temperature limit	Up to 100°C without sleeve
Pathlength	5 mm
Probe material	Stainless steel outer sleeve, aluminum second surface mirror
Pressure	100 psi
Fiber jacketing	Black PVC monocoil
Connectors	SMA 905

### PEEK pH Probe

Engineering Specifications	TRANS-PH-KIT-PEEK
Probe diameter	6.35 mm (1/4")
Probe length	107.9 mm
Fiber length	2 meters
Breakout	1.5 meters from the end of the probe
Temperature limit	Up to 200°C with PEEK sleeve
Pathlength	Adjustable (2 – 10 mm)
Probe material	PEEK
Fiber jacketing	Black PVC monocoil
Connectors	SMA 905

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## Thin Film Metrology

### Spectroscopic Reflectometry Systems

#### Sample NanoCalc Applications

- \* Transmission and reflection measurements of anti-reflective and hardness coatings
- \* Analysis of medical coatings and catheter balloon foils
- \* Testing of the hardness and wear of coatings
- \* Measurement of the thickness of thinned silicon wafers
- \* Determination of photoresist layers for masks
- \* Analysis of coatings applied for weather or dirt resistance (Lotus Effect)
- \* Measurement of coatings inside beverage containers
- \* Air gap measurements
- \* Analysis of optical disk coatings

#### NanoCalc Features

- \* Resolution to 0.1 nm
- \* Able to measure stacks of up to 10 layers
- \* Thickness and refractive index data
- \* Sophisticated algorithms for defect and roughness tolerance measurements
- \* Large database to ensure accuracy of a broad range of materials
- \* Adapters for complex geometries and accessories for thickness mapping

#### Specifications

Item	NANOCALC-VIS-PRECON	NANOCALC-XR-PRECON	NANOCALC-DUV	NANOCALC-NIR
Wavelength range	400 - 850 nm	250 - 1050 nm	~200 - 1100 nm	900 - 1700 nm
Thickness range	50 nm - 20 μm	10 nm - 100 μm	1 nm - 100 μm	100 nm - 250 μm
Optical resolution	0.1 nm	0.1 nm	0.1 nm	0.1 nm
Repeatability	0.3 nm	0.3 nm	0.3 nm	1.0 nm
Angle of incidence	90°	90°	90°	90°
Number of layers	Up to 10	Up to 10	Up to 10	Up to 10
Refractive index	Yes	Yes	Yes	Yes
Test materials	Transparent or semi-transparent thin film materials	Transparent or semi-transparent thin film materials	Transparent or semi-transparent thin film materials	Transparent or semi-transparent thin film materials
Reference needed	Yes (bare substrate)	Yes (bare substrate)	Yes (bare substrate)	Yes (bare substrate)
Measurement modes	Reflection and Transmission	Reflection and Transmission	Reflection and Transmission	Reflection and Transmission
Rough materials capable	Yes	Yes	Yes	Yes
Measurement speed	100 ms to 1 s	100 ms to 1 s	100 ms to 1 s	100 ms to 1 s
On-line capable	Yes	Yes	Yes	Yes
Height adjustment	with COL-UV-6.35 (10-50 mm)	with COL-UV-6.35 (10-50 mm)	with COL-UV-6.35 (10-50 mm)	with COL-UV-6.35 (10-50 mm)
Spot size	200 μm or 400 μm standard; 100 μm available upon request	200 μm or 400 μm standard; 100 μm available upon request	400 μm standard; 200 μm available upon request	400 μm standard; 200 μm available upon request
Microspot	Yes (w/microscope)	Yes (w/microscope)	Yes (w/microscope)	Yes (w/microscope)
CCD color	Yes (w/microscope)	Yes (w/microscope)	Yes (w/microscope)	Yes (w/microscope)
Mapping option	150 mm (6") and 300 mm (12") xy-scanning stages	150 mm (6") and 300 mm (12") xy-scanning stages	150 mm (6") and 300 mm (12") xy-scanning stages	150 mm (6") and 300 mm (12") xy-scanning stages
Vacuum capable	Yes	Yes	Yes	Yes

Item	Wavelength Range	Optical Layer Thicknesses
NANOCALC-VIS-PRECON	400 - 850 nm	50 nm - 20μm
NANOCALC-XR-PRECON	250 - 1050 nm	10 nm - 100μm
NANOCALC-DUV	~200 - 1100 nm	1 nm - 100μm
NANOCALC-NIR	900 - 1700 nm	100 nm - 250μm

#### Specifications Standard Operating Software (Required)

Item	Description	Required for
NANOCALC-10-N	Thin film measurement software for Windows; measurement, simulation and analysis of up to 10 layers; refractive index analysis also possible (refractive index for multiple layers requires SCOUT software)	Any NanoCalc measurement of thin films up to 10 layers

## Add-on Software Options

Item	Description	Required for
SCOUT-FULL VERSION	Spectrum simulation program for Windows XP/7 (32 bit). Computes reflectance, transmittance, absorbance or ellipsometry spectra and fits your model to measured spectra by manual, graphical or automatic parameter variation. SCOUT can be controlled by OLE automation controllers.	Spectral simulation; use w/NANOCALC-10-N
NANOCALC-MAPPING	Mapping module software (needs NANOCALC -10-N) is 3D-mapping module with control of 150 mm and 300 mm mapping stages	Systems using mapping stages
NANOCALC-ONLINE	On-line module software (needs NANOCALC-10-N) includes external trigger and manual multipoint measurement with data transfer into Excel, plus statistical feature and 1D-Plot; also provides online display of XY-graphs and histograms of layer thickness and removal rate	On-line applications
NANOCALC-MULTIPOINT	Multipoint module software (needs NANOCALC -10-N) for manual multipoint measurement; provides result list with data transfer into Excel, CSV data, statistical and 1D-graphic; lets you measure on a mouse click, keyboard key or external trigger	Multipoint measurements
NANOCALC-REMOTE	Remote module (needs NANOCALC -10-N). Active-X functionality allows control of most NanoCalc functions from any other software.	Controlling NanoCalc functions w/other software

## Fiber Assemblies for NanoCalc Systems

Item	Description	Use with NanoCalc Models
NC-2UV-VIS400-2	Bifurcated UV-VIS fiber, 400 $\mu\text{m}$ , 2 m long, 2 x SMA 905 connectors, flexible metal jacketing, common end with stainless steel ferrule 6.35 mm x 50 mm	NANOCALC-VIS NANOCALC-XR
NC-7UVS400-2	NANOCALC-DUV reflection probe, 400 $\mu\text{m}$ , 2 m long, 6 illuminated fibers, flexible metal jacketing, common end with stainless steel ferrule 6.35 mm x 50 mm	NANOCALC-DUV
NC-7VIS-NIR400-2	NANOCALC-NIR reflection probe, 400 $\mu\text{m}$ , 2 m long, 6 illuminated fibers, flexible metal jacketing, common end with stainless steel ferrule 6.35 mm x 50 mm	NANOCALC-NIR

## Reflection Probes for NanoCalc Systems

Item	Description	Use with NanoCalc Models
NC-7UV-VIS200-2	Reflection probe, 6 illumination, 1 read fiber, 200 $\mu\text{m}$ UV-VIS fibers, 2 m long, flexible metal jacketing, stainless steel ferrule 6.35 mm x 50 mm, 2 x SMA 905 connectors	NANOCALC-VIS NANOCALC-XR
NC-7UV-VIS200-2-SMA	Reflection probe for use w/MFA-PT microscope adapter, 6 illumination, 1 read fiber, 200 $\mu\text{m}$ UVVIS fibers, 2 m long, with flexible metal jacketing, stainless steel ferrule 6.35 mm x 50 mm, 3 x SMA 905 connectors	NANOCALC-VIS NANOCALC-XR

## Accessories

### Reference Wafers

We offer two Si-SiO<sub>2</sub> reference wafer options for measuring the thickness of thin, transparent films on substrates such as silicon wafers and optical layers. The STEP-WAFER covers UV-VIS wavelengths and the STEP-WAFER-600-1100 covers VIS-NIR wavelengths. Each wafer is a 100mm diameter, 5-step wafer with calibrated thickness range of 0-500nm or 600-1100nm.

## Thin Film Metrology

### SpecEI-2000 Ellipsometry System

The SpecEI-2000 is a benchtop thin film measurement system utilizing spectroscopic ellipsometry to measure multilayer, semi-transparent samples such as flat wafers or glass plates. The Spec-EI-2000 is affordable, compact (52 cm x 33 cm x 24 cm) and convenient, with easy placement of the sample and one-button operation. SpecEI is part of the Mikropack line of thin film metrology systems. Models are available for 300-1000 nm (SPECCEL-2000-UV-VIS-NIR) and 400-1000 nm (SPECCEL-2000-VIS-NIR). SpecEI systems include an integrated spectrometer, broadband light source and controller; software is available separately.

#### SpecEI Features

- \* Film thickness accuracy to 1.0 nm
- \* Resolution to 0.1 nm
- \* Multi-layer stack measurements up to 25 layers
- \* Single film thickness up to 10  $\mu\text{m}$
- \* Spectral ranges from 300-1000 nm
- \* Standard spot size 0.4 mm x 1.2 mm
- \* Ideal for flat, semi-transparent samples such as wafers, glass, films and foils
- \* 3D mapping, reference wafers, accessories and other options available
- \* Accompanying software allows generation and recall of measurement recipes for one-step, repetitive measurements
- \* Accessories for thickness mapping

#### Software for SpecEI-2000 Systems

Powerful SpecEI software offers a range of modeling possibilities such as Cauchy, OJL, Tauc-Lorentz, Drude, EMA and different types of oscillators. The software also stores specific measurement routines, reducing the tedium of repetitive measurements and easing integration.

#### Specifications

System Performance	
Thickness range	1 nm - 10 $\mu\text{m}$
Resolution	0.1 nm
n and k analyzer	Values calculated for complete spectral range
Mathematical models	Extensive range of options includes constant refractive index, harmonic oscillator and imported dielectric functions
Measurement speed	7 - 13 seconds
Repeatability	70 nm for $\text{SiO}_2$ on Si, $\cos(\Delta) \pm 0.0003$ , $\tan(\Psi) \pm 0.0002$
Sample size	Desktop up to 150 mm diameter; mapping up to 300 mm diameter
Sample thickness	5 mm (maximum)
Optical	
Wavelength range	(UV-NIR) 300 - 1000 nm or (VIS-NIR) 400 - 1000 nm
Optical resolution	1.0 nm
Beam diameter	400 - 1200 $\mu\text{m}$
Angle	70°
Computer	
Software	Windows XP/7 (32bit) software; also, recipe structure, administrator-user compatible
Hardware	PC with Windows XP/7 (32bit)
Standard Operating Software (required)	
SPECEL-ELLICALC	User-friendly 32-bit Windows software for ex situ direct measurement of thickness and n and k coefficients; recipe structure and administrator/user capabilities
SCOUT-FULL VERSION	SCOUT software. Spectrum simulation program for Windows XP/7 (32 bit). Computes reflectance, transmittance, absorbance or ellipsometry spectra and fits your model to measured spectra by manual, graphical or automatic parameter variation. SCOUT can be controlled by OLE automation controllers.
Add-on Software Option	
SPECEL-MAPPING	Mapping module software (must be used with SPECEL-ELLICALC software). Complete mapping module with software control of 150 mm and 300 mm mapping stages; includes XYZ controls and 3D-graphics

#### Accessories for SpecEI-2000 Systems

We offer reference wafers and mapping stages for use with SpecEI systems. Mapping stages must be purchased at time of SpecEI system order. Replacement parts are also available.

Item	Description
MAPPING-6-INCH-SE	150 mm x 150 mm XY-scanning stage; fully automatic w/vacuum chuck, 2 motors with encoders, control system integrated 2-axis CNC controller, RS-232 interface; portal structure
MAPPING-12-INCH-SE	300 mm x 300 mm XYZ-scanning stage; fully automatic w/vacuum chuck, 2 motors with encoders, control system integrated 2-axis CNC controller, RS-232 interface; portal structure
STEP-WAFER	Reference Si-SiO <sub>2</sub> -step-wafer, 5 steps from 0-500 nm, calibrated, 100 mm diameter
STEP-WAFER-600-1100	Reference Si-SiO <sub>2</sub> -step-wafer, 5 steps from 600-1100 nm, calibrated, 100 mm diameter

## Thin Film Metrology

### PlasCalc-2000 Plasma Monitoring and Control System

The PlasCalc-2000 system provides real-time, in situ analysis of the optical emission spectra acquired during plasma processes. The system has all the tools necessary for monitoring and controlling a running process, with sophisticated algorithms for data acquisition and signal treatment. Wavelength, recipe and formula editors allow quick creation of efficient recipe functions for data handling and easy combination of many recipes to build comprehensive strategies for system response during monitoring. PlasCalc-2000 (item PLASCALC-2000-UV-VIS-NIR) is part of the Mikropack line of thin film metrology systems.

The PlasCalc-2000 is 257 mm x 152 mm x 263 mm and includes a spectrometer (200-1100 nm) with 16-bit D/A converter and 12 VDC power supply. Operating the system requires PlasCalc software and 400  $\mu$ m optical fibers, which are available separately.

### Sample Plasma Monitoring Applications

- \* Plasma etching
- \* Plasma chamber health control
- \* Planarization of blanket polysilicon
- \* Application of protection coatings
- \* Abnormal process phenomena
- \* Process optimization
- \* Surface cleaning processes

### Spectroscopic

Spectral range	200 - 1100 nm
Optical resolution	1.0 nm (FWHM)
Fiber optic connector	SMA 905

### Electronics

Power supply	12 VDC, 1.25 A
Power requirements	90 - 240 VAC, 50/60 Hz
D/A converter	14 bit " I/O "
Digital I/O	8x TTL
Analog output / voltage sign	4x [0-10V]

### Computer

Software	PlasCalc software (basic operating software); also, SpecLine atomic emission line analysis software (add-on package)
Hardware	PC with Windows XP/7 (32 bit)

### Software for PlasCalc-2000-UV-VIS-NIR System

#### Standard Operating Software (Required)

Item	Description
PLASCALC-SOFTWARE	Real-time and in situ access to all functions via easy menu-oriented software interface. Multiple plasma species can be picked by a mouseclick, mathematically calculated, subtracted and monitored in real time and in situ. Capabilities: four digital inputs/outputs, 4 analog output channels; display messages with tone signal can be set up for end point, start, stop and limit.

#### Add-on Software Option

Item	Description
SPECLINE	Software for spectroscopy, astrophysics, plasma science or plasma processing. This tool supports and makes it easy to evaluate spectral data, e.g. finding specific lines in spectra, identifying unknown peaks, identifying atomic lines and molecular bands or comparing data from different measurements in spectra data. Almost instantly detects peaks of lines and bands using several powerful filter functions. Extensive database for atoms and molecules included; only for Windows XP/7 (32 bit).

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## Software

### SpecLine Software for Compound ID

Flame Series, USB Series, HR Series, QE Pro Series, and MAYA2000 Pro Series.

#### Identifies Elements and Compounds

SpecLine Software is a powerful tool designed for identifying atomic emission lines and molecular bands in spectral data. SpecLine's advanced evaluation, search, compare and identify functions and its extensive library of over 100 elements and over 400 compounds enable you to quickly identify unknown lines, peaks and bands. SpecLine was designed for scientists, engineers and researchers using emission spectroscopy in fields such as astrophysics, the plasma sciences and plasma processing.

#### Searching and Comparing Data

In the Line Identification window, you can define all the parameters for your search in a Periodic Table screen (top right), and begin the process with just a single click. SpecLine can analyze even the most complex spectral data, including spectra with double lines, line shoulders and complex band structures. Up to 12 separate spectra, even if they are in different file formats, can be combined for comparative purposes.

#### Identification


SpecLine applies a variety of sophisticated filter functions such as Wiener-Fourier and polynomial noise removal to identify the elements and compounds in your spectra. After SpecLine applies comparative searches to its extensive atomic, ionic and molecular database of over 100 elements in several ionization states and over 400 elemental compounds, it provides detailed data on each identified peak and line, such as the name of the element, the peak's wavelength, the electron voltage and its transition state and quantum number.

#### Opening Spectra and Saving Data


SpecLine can analyze spectral data from various spectroscopy software applications; it also can open all Ocean Optics software file formats as well as SPC and ASCII file formats. In addition, you can save all of SpecLine's identification data in its native file format or export it into various applications, such as Excel.

#### Hardware Keys


SpecLine Software comes with a USB or parallel-port (printer port) hardware key. The key is a security device to protect against unlicensed copies. It connects to an input/output port on your computer and must be used to run the software. SpecLine-U comes with a USB dongle for use with Windows XP, Vista and Windows 7 (32-bit). SpecLine-P comes with a feed-through parallel port hardware key.

In the Identify Lines window, you can search atoms from a Periodic Table, molecules from an extensive list, and elements in single or multiple ionization states.

In this window, a search on the peaks and lines in a spectrum has been completed and identified successfully.

This window demonstrates SpecLine's ability to provide detailed data on just one emission line.

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